

# HOME OF RHEOLOGY



GÖTTERT  
 THIS IS RHEOLOGY

# MELT FLOW INDEXER

Modular Melt Indexer series and options

Technical modifications / Highlights new MI series compared to previous models

mi1 | mi2 | mi3 | mi40 | MI-ROBO

Launched at the K2025 trade fair - Oct25

# MELT FLOW INDEXER

ISO 1133 | ASTM D1238 | D3364



mi1




mi2




mi3



mi40



MI-ROBO



# MELT FLOW INDEXER **mi1**

Highest precision – minimal automation!

- **NEW – Optionally selectable WSVGA Touchscreen 4,33"**
- **NEW** - Standardization of test chamber – easy change of test barrel for different applications (corrosive/abrasive) – **Proven Technology, as already used with the models mi3 and mi40 models**



# MELT FLOW INDEXER **mi2**

Highest precision – half automated – equipped with test-weight Lifting system – single test weight measurements



- **NEW** - WSVGA Touchscreen 7"
- **NEW** - Displacement sensor 0,003mm/ impuls - **previously 0.025mm/ impuls**
- **NEW** - Standardization of test chamber – easy change of test barrel for different applications (corrosive/abrasive) – **Proven Technology, as already used with the models mi3 and mi40**
- **NEW** - Automatic pre-positioning of the piston by the new lifting system and parameter settings - **No material sensor required, as with the previous model 2.2**



# MELT FLOW INDEXER **mi3**

Highest precision- incl. test-weight magazine - manual test weight selection - single test weight measurements



- **NEW** - WSVGA Touchscreen 7"
- **NEW** - **optionally selectable** Force controlled **compressing /ejection up to 80 kg**. The test-weight lifting motor can be used (position and time controlled) to pre-load the sample before the test or eject the material after the test. The maximum available force is 80 kg as a combination of the 21,6 kg load and the motor at full force. - **Not possible with the previous model** -

**NEW** - **optionally selectable** Semi-Automatic Cleaning Unit

For barrel cleaning without compressed air. The module is firmly mounted on the test device. In addition to the up/down movement the rotating cleaning piston guarantees a reliable and efficient cleaning of the test barrel. - **Not possible with the previous model** -

# MELT FLOW INDEXER **mi40**


Highest precision and automation - incl. test-weight magazine - automatic test weight selection - single and multi test weight measurements - viscosity option



- **NEW** - WSVGA Touchscreen 7"
- **NEW** - Force controlled **compressing/ejection up to 80 kg**. The test-weight lifting motor can be used (position and time controlled) to pre-load the sample before the test or eject the material after the test. The maximum available force is 80 kg as a combination of the 21,6 kg load and the motor at full force.  
**With the previous model, only possible up to 60 kg**

# MELT FLOW INDEXER **MI-ROBO**

Highest precision - Fully automated Melt Flow Indexer



- **NEW** - Determination of Melt Density (Method A)

The MI-ROBO optionally enables the automatic determination of melt density according to Method A in compliance with ISO 1133 / ASTM D1238 Procedure A

In this process, the melt volume is recorded based on the melt index, and the material density is automatically calculated. This method increases reproducibility and minimizes operator errors.

- **NEW - now available** Visualization and Control software miCONNECT



# GÖTTFERT



THIS IS RHEOLOGY

Martin Ackermann - 28.10.2025



Connect with us on  
**LINKEDIN**



Follow us on  
**TWITTER**



Watch us on  
**YOUTUBE**



Like us on  
**FACEBOOK**



Find us on  
**INSTAGRAM**